

## eAirway.04 - Airway Device Placement Confirmed Method

## Definition

The method used to confirm the airway device placement.

National Element	No	Pertinent Negatives (PN)	No
State Element	Yes	NOT Values	Yes
Version 2 Element		Is Nillable	Yes
Usage	Recommended	Recurrence	0 : M

## Associated Performance Measure Initiatives

Airway

## Attributes

**NOT Values (NV)**

7701001 - Not Applicable

7701003 - Not Recorded

7701005 - Not Reporting

**CorrelationID**

**Data Type:** string

**minLength:** 2

**maxLength:** 255

## Code List

Code	Description
4004001	Auscultation
4004003	Bulb/Syringe Aspiration
4004005	Colorimetric ETCO <sub>2</sub>
4004007	Condensation in Tube
4004009	Digital (Numeric) ETCO <sub>2</sub>
4004011	Direct Re-Visualization of Tube in Place
4004013	Endotracheal Tube Whistle (BAAM, etc.)
4004015	Other
4004017	Visualization of Vocal Cords
4004019	Waveform ETCO <sub>2</sub>
4004021	Chest Rise

## Data Element Comment

If the invasive airway is confirmed via multiple methods each method should be documented individually by time, method, and type of individual.

## Version 3 Changes Implemented

Added to better document airway management.

## Associated Validation Rules

Rule ID	Level	Message
nemSch_e001	Error	When Airway Device Placement Confirmed Method is empty, it should have a Not Value (Not Applicable, Not Recorded, or Not Reporting, if allowed for the element) or a Pertinent Negative (if allowed for the element), or it should be omitted (if the element is optional).
nemSch_e002	Error	When Airway Device Placement Confirmed Method has a Not Value (Not Applicable, Not Recorded, or Not Reporting), it should be empty.
nemSch_e009	Warning	When Airway Device Placement Confirmed Method has a Not Value, no other value should be recorded.